





An effective, chlorine-free sanitiser for baby's bottles, teats and associated utensils.

Unlike chlorine based products Little One Sanitiser will not form harmful by-products, as it decomposes to water and oxygen only.

- The 'must-have' item of modern parents. It can be used for an extraordinary range of tasks including:
  - sanitising your child's bottles, teats and utensils.
  - children's toy sanitising.
  - sanitising drink bottles / sports bottles.
  - sanitising cosmetic tools such as hair combs, tweezers etc.
  - freshening cut flowers.
- A truly versatile and innovative improvement on market equivalents.

FEATURES	BENEFITS
Chlorine-free formula	<ol> <li>Will not form harmful by-products</li> <li>Active hydrogen peroxide decomposes to water and oxygen only</li> <li>Will not produce obnoxious tastes or odours</li> <li>Safe to use in septic or multi-stage treatment systems</li> </ol>
Versatile	<ol> <li>Can be used to sanitise and to freshen cut flowers</li> <li>Sanitise cosmetic tools such as tweezers</li> <li>Multi-purpose - a great advantage over tablet-style equivalents</li> </ol>
Economical	Price-per-use compares favourably with popular supermarket brands     Outperforms most tablet-style products in price-per-use
Safe	Price-Harmless on all surface - will not bleach or damage

## **DIRECTIONS FOR USE**

*Sanitising baby bottles, teats, utensils, drink bottles etc:* Dilute 10mL per 1Lt water. Soak for a minimum of 60min. Solution remains active for up to 24 hours.

*Children's toys:* (To soak) Dilute 10mL per 1Lt water. Soak for a minimum of 60 mins. Ensure that toys are thoroughly dried (inside) to eliminate the potential for bacteria to grow inside a toy.

(To wipe) Dilute 10mL per 1Lt water. Use clean cloth/sponge to dip in solution and wipe over exterior of toys. Ensure that the toy is allowed to dry thoroughly.

*Sanitising cosmetic tools such as hair combs, tweezers etc:* Dilute 100mL with 500mL water. Soak for a minimum of 60mins. Solution remains active for up to 24 hours.

Freshening cut flowers: Add 2mL to a vase of water. Change solution every 3-4 days.

**INGREDIENTS** 

Activated Hydrogen Peroxide 6%

